Gas fermentation of *Cupriavidus necator* variants at elevated pressure

### JASBIR SINGH and GERARD GARDNER (HEL LTD), TUCK SENG WONG and KANG LAN TEE (University of Sheffield), STEPHANE E GUILLOUET (LISBP, INSA, France)





better chemistry - faster

C1net Conference 21 -23 January, 2019





## **Gas Fermentation Rate – commercial viability**

**<u>C1-based</u>** fuels and more valuable organics need to be produced at rates which are fast enough to justify the investment in plant and people: presently big challenge for large scale production

### **Reasons**

- Low density of gases (~ 1/1000<sup>th</sup> of liquids)
- Low gas solubility (CH4,  $H_2$  and  $O_2 \sim 1$  mmole/kg water,  $CO_2 \sim 30$  mmole/kg water)

C1 gas transfer rate is slow and limits conversion rate.

## **Solution**

Speed up gas transfer



## **Gas Concentration in Broth**

Mass balance around bioreactor describes changes in dissolved oxygen concentration C (in milligrams per liter) in the nutrient broth gives:





# **Gas Transfer into Liquid**

### <u>Agitation and Sparging – first step</u>

The Mechanical design of agitator blades, speed of agitation and break up into small bubbles provides an important step in getting the gas to the bacteria/cells: expressed in terms of kLa

### <u>Solubility</u> – second step

Molecular forces between the gas and liquid, expressed in terms of solubility, then determines how much gas is actually retained in the liquid.

Hence: Gas uptake rate =  $La (C^* - C)$ 

## How to increase solubility C\* ?





Alternate oxygen and nitrogen purging at different pressure. Oxygen saturation level increases with pressure.



### Maximum Gas transfer rate into Liquid

### Agitation speed, impeller and sparging fixed



Gas transfer rate increases linearly with pressure



### Typical Systems used for C1 Fermentation





## Carbon utilization of *Cupriavidus necator*

Microorganism widely recognized for conversion of C1 sources into a range of different organics, according to the choice of bacterial strain and features of bioengineering.

Three examples:

- C. Necator Re2133 (studies of strain growth with CO2/H2/O2 feed)
- C. Necator pEG7d (studies of autotropic isopropanol production)
- C. Necator H16 (studies of strain growth with both gluconate and CO2 as carbon sources)



#### Autotrophic growth of STRAIN in bioreactor

#### Autotrophic growth of Cupriavidus necator Re2133 at atmospheric pressure

Grown in stirred reactor by supply of CO2/H2/O2

Fermentation stopped when the maximum O2 transfer capacity is reached (DO near zero)

maximum O2 transfer rate determines maximum growth rate





#### Autotrophic production of strain – at pressure

Autotrophic growth of Cupriavidus necator Re2133 at high pressures

Efforts to increase growth rate by elevated pressure of gases

Pressure ramps to continue growth – up to 4 bar

Growth rate does not increase but growth is sustained for longer OD now ~ 1.2 instead of ~ 0.2 at atmospheric pressure



### **Autotrophic production of isopropanol**





#### Autotrophic production of isopropanol





### **Isopropanol Production - benefits of elevated pressure**

#### > Maximum rate of IPA production in CO<sub>2</sub> matched rate with fructose

Production phase on	[isop] <sub>max</sub> (g L <sup>-1</sup> )	qIsop <sub>max</sub> (Cmol Cmol <sup>-1</sup> h <sup>-1</sup> )	Specific production (g isop/g <sub>CDW</sub> )
<b>CO<sub>2</sub></b> from Gas Mix	0.25	0,035	0.25
CO <sub>2</sub> from Ind. Gas	2.2 (3)	0.07	1.0
Fructose	9.0 (16)	0,062	2



### Fermentation of *Cupriavidus necator* H16

### **Objective**

#### C. necator H16 during fermentation:

(a) under heterotrophic cultivation conditions(b) under autotrophic cultivation conditions

#### Typical gas feed composition:

 $85 \% H_{2,} 5 \% O_{2'} 10 \% CO_2$ 



#### Heterotrophic growth of *C. necator* H16 at 1 bar in 1% sodium gluconate

#### Reference Run at 1bar



Heterotrophic growth of *C. necator* H16 at 1 bar in 1% sodium gluconate. The maximum OD<sub>600</sub> (=15) and calculated  $\mu_{max}$  (=0.326 h<sup>-1</sup>) are similar to values achieved in tube cultures.



#### Heterotrophic and Autotrophic cultivation of C. necator H16



Growth curve of C. necator H16 in gas fermentation



### Fermentation at 1-4 bar pressure





### Cupriavidus Necator pEG7d at elevated pressure

Cell's viability determination by Flow Cytometry



Increase in Pressure has no negative impact on cell's viability (within the range 0 to 5.3 bars overpressure)



## CONCLUSIONS

**<u>1.</u>** Limits of gas solubility and density can be overcome by operation of reactor at elevated pressure.

**<u>2</u>**. Bioreactor systems operating up to 10bar have been developed and successfully operated.

<u>3.</u> Specific Bacterial Strains need to developed to survive high pressure conditions.

<u>4</u> Cupriavidus Necator pEG7d viable at least to 5bar (GUILLOUET, LISBP, INSA, France). Already able to produce alcohols from C1 gas to over 2g/L

<u>4</u> Cupriavidus Necator H16 at present ~ 3bar, under development (Wong and Tee at Sheffield Univ., UK).

